



The Oxygen Recovery Mask is indicated whenever there is a need to deliver supplemental oxygen therapy to an animal such as a dog or a cat.

A. Procedure

Animal Care

1. Prepare equipment:
 - a) Oxygen Recovery Mask (select appropriate size)
 - b) Oxygen tubing
 - c) Connect the oxygen supply line to the Oxygen Recovery Mask.
 - d) Oxygen fitting maybe removed to allow connection of a BVM.
2. Seal mask over animal's mouth and nose.
3. Set appropriate oxygen flow rate to at least 10 L/min.
4. Monitor animal to assure animal does not run off with mask attached.

B. Special Considerations

1. Masks come in three sizes: small, large, and feline.
2. Masks are reusable. Wash masks with soap and water. Air dry.
3. If the animal is combative do not force the mask on the animal.