



CHEATSHEET

This Controlled Substance eLog (CS eLog) Cheatsheet is a guide governed by and is subordinate to Procedure 37 and at no time supersedes its authority or its direction.

3 Functions of the CS eLog

OIC Exchange: The transfer of the CS eLog ownership at morning shift change or at any OIC exchange of responsibility.

Medication Administration: Logging patient controlled substance administrations.

Unit Transfer: Transferring controlled substances to another MDRF unit.

Things to Know !

Access CS eLog from Tablet?

On the Tablet → CS eLog Icon

Access CS eLog from Station Computer

MDFRNET → Controlled Drugs eLog

Login & Witness Signatures?

MDC: (User name / password)

OIC Exchange


When to do an OIC Exchange?

- Every morning @ face-to-face, exchange of responsibility.
- @ Every shift change including (holdover, OT) when on the same unit.
- ASAP when off-going, face-to-face is not possible.

Who logs in to the CS eLog at the OIC exchange?

- The On-coming OIC.

What to do?



- Click,  Exchange
- Count all available narcotics.
- Enter counts in the CS eLog.
- Get Witness signature.

Who Witnesses / Signs?

- Off-Going OIC
- (or) Units' Paramedic, only if OIC face-to-face is not possible.

Medication Administration

How to enter?

1. Click,  Administer
2. Click,  Add – next to the pre-populated ePCR medication entry, to auto-fill fields. (if available)
3. Verify/Correct info. entered.(edit if req.)
4. Get Paramedic / Witness signature and click save.

Medication given & Not Pre-populated?

1. Contact an EMS Captain to do an Override
2. Add "Narrative Addendum" to EMS report if needed.
3. EMS Captain will make the entry in the units log as an Override.

Unit Transfer

 Transfer

When to do a transfer?

- Any time a MDRF unit transfers controlled medications to another unit.

Who can do a transfer?

- Any MDRF unit.

Who Logs the transfer and Where?

- The sending MDRF unit documents the transfer in their log.

Who Witnesses / Signs?

- The receiving units OIC must witness the sending units CS eLog.